sartorius stedim
biotech

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Internet: www.sartorius-stedim.com
E-mail: tech.pubs@sartorius-stedim.com

Printed in Germany on chlorine-free paper.

If you require mandatory information regarding specific properties of your reverse osmosis system that goes beyond that which is described here, please contact us. Reprinting, translation and copying in any form, even in excerpts, requires the written permission of Sartorius Stedim Biotech GmbH.

arium® is a registered trademark of Sartorius Stedim Biotech GmbH

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Table of Contents

1. General	4
1.1 Notes Regarding this Manual	4
1.2 Safety Information	5
1.3 Terms	5
1.4 Intended Use	5
1.5 Product Description	6
1.6 Operating Data	7
1.7 Technical Specifications	7
2. Initial System Operation	8
2.1 System Placement	8
2.2 Water Quality	8
2.3 Water Connections	9
2.4 Electrical Connection	9
2.5 System Startup	10
2.6 Adjustment of Concentrate Flow Rate	11
2.7 Setting the Permeate Flow Rate	12
3. System Operation	13
3.1 Automatic	13
3.2 Maintenance	13
3.3 Trouble-Shooting	13
3.3.1 System Capacity has Dropped	13
3.3.2 The Quality of the Treated Water is Higher than 20 µmS/cm	14
3.3.3 Alarm: Lack of Input Pressure	14
3.3.4 Alarm: Transport Pump	15
3.3.5 Alarm: Tank Fill Level Too Low	15
3.3.6 Alarm: Tank Fill Level Too High	15
3.3.7 System Will Not Start	15
3.3.8 Transport Pump Turns On Off	16
3.3.9 Alarm Display	16
3.3.10 Time Setting	16
3.4 Cleaning Replacing Membranes	17
3.4.1 Cleaning Procedure	17
3.4.2 Tank Cleaning Procedure	18
3.4.3 Replacing Membranes	18
4. Appendix	21
4.1 Flow Diagram	21
4.2 Detailed Systems Drawings	22
4.2.1 Layout Overview	22
4.2.2 Front Side View	23
4.2.3 Back Side View	24
4.2.4 Top View	25
4.3 Wiring Diagram	26
4.4 Main Board	27
4.5 Service and Maintenance of RO System Maintenance Schedule	29
4.6 Operating Log	31
4.7 Disposal	33
Contacts	67

1. General

1.1 Notes Regarding this Manual

- ▶ Please read this entire manual carefully before using the device.
- ▶ Please ensure that you read the safety information carefully.
- ▶ This manual is considered a part of the product. Keep it in a safe and easily accessible location.
- ▶ You can request a replacement copy if it is lost.

Symbols and Signs

The following symbols are used in this manual:

This indicates danger that is immediate or can always be expected which would lead to serious injury or death.

This indicates danger that has a high probability of occurring which would lead to serious injury or death.

This indicates danger that can lead to medium and/or light injury.

This indicates danger that may cause property damage.

This symbol indicates useful information and tips.

- ▶ Indicates a required action
- ▷ Describes the result of an action
- 1. If a procedure has multiple steps...
- 2. ... the steps are numbered consecutively.
- Indicates an item in a list

1.2 Safety Information

The RO 613L reverse osmosis system has been constructed in accordance with the European Directives as well as international regulations and standards for operation of electrical equipment, electromagnetic compatibility and applicable safety requirements. Improper use or handling, however, can result in damage and/or injury. The equipment may only be operated by trained personnel. Operators must read these installation instructions, particularly the safety information, and must be familiar with the operation of the equipment. The operator must supplement these warnings and safety precautions if necessary and brief operating personnel accordingly.

Before operating the RO 613L reverse osmosis system, be sure to read the "Safety Information" section in this operating manual.

1.3 Terms

The following technical terms are used in this operating manual:

Softening system	An optional system that removes calcium carbonate and magnesium from hard water.
Concentrate	The water that is led to the outlet. This water contains the salts and minerals that have been removed from the raw water.
Conductivity	The amount of electrical conductivity of the water being treated (permeate), measured in $\mu\text{mS/cm}$. The lower the value the better the water quality.
Level sensor	A float that turns the transport pump on/off and stops the pump, if necessary, to prevent the storage reservoir from emptying completely.
RO membranes	System filters that remove salt from the raw water under high pressure and flow.
Permeate	The treated, totally desalinated water produced by the RO system and sent to the tank.
Reservoir pump Transport pump	The integrated pump that supplies treated water from the system tank to the user.
RO	This is an abbreviation for reverse osmosis.
Raw water	The water supplied to the RO system to be desalinated.
TDS	Totally dissolved salts, measured in mg/l .

1.4 Intended Use

- The following points must be followed to ensure proper installation and operation of your RO system:
1. The included startup log should be filled out during startup and archived together along with the operating log.
 2. The operating log (see appendix) must be updated on a regular basis.
 3. A non-pressurized floor drain must be located near the system.
 4. The RO system removes 95-98% of all salt from the raw water. Secondary treatment using a mix-bed or similar, for example, may be necessary if better water quality is required.
 5. This operating manual was developed according to regulations and legal guidelines currently available at the time.
It is structured so that it can be implemented and used by trained personnel. Activities that go beyond the scope of that which is described here are only to be carried out by SSB service personnel.
 6. The warranty becomes null and void if the initial startup of the system is not carried out by an authorized service technician.

1.5 Product Description

The arium RO613L system was designed to create reverse osmosis water.

The water supplied to the system must be drinking water quality and may not contain any residue chlorine. If there is more than 0.1 mg/l of free chlorine in the raw water, you must add an activated carbon filter to the system (accessory).

The flow rates and qualities listed are based on drinking water quality with a raw water salinity of a max. of 1000 mg/l TDS and 10°C. If you have raw water quality that deviates, we ask that you please contact your supplier. See also section 2.2. Water Quality.

The raw water is pumped through the RO membrane(s) via the high pressure pump whereby the salts contained in the water are retained and only permeate passes through the fine openings of the membranes and on to the tank. The water with the concentrated salts (concentrate) is guided through the RO module(s) directly to the drain. The ratio between permeate and drainage is regulated manually via needle valves.

The raw water passes through a 5 µm prefilter cartridge. If required, a filter with an activated carbon cartridge (chlorine content > 0.1 mg/l) must be installed.

The pressure switch on the supply water input is used to control the water input pressure. The water input pressure must be > 0.5 bar. If the pressure is too low, the high pressure pump turns off.

The tank is equipped with a level switch that regulates the raw water intake and the operation of the high pressure pump. In addition, the transport pump is also turned off before the tank is emptied completely.

The permeate created is fed to the point of consumption automatically via the transport pump. The maximum flow rate is limited by the permeate capacity of the system and the size of the storage tank (100 l).

Optional accessories can also be used, if required, to increase the water quality, e.g. an ion exchanger (mixed-bed) and/or a UV lamp (for sterilization).

Your RO system is mounted on rollers and, therefore, should be attached to flexible connections. The system can be placed under a table, for example, and still remain easily accessible for use.

The RO membranes have a long service life under normal operating conditions. However, even good raw water quality may cause a film of impurities to build up that can reduce the permeate yield little by little. When the yield has been reduced by 10%, the membranes must be cleaned. Regular cleaning and observing the prescribed maintenance schedule will maintain original system performance over the long-term.

1.6 Operating Data

RO 613L SYSTEMS	613L050	613L100	613L150	613L200	613L250	613L300
Capacity [l/h]**	50	100	150	200	250	300
Max. water recovery [%]***	40-80	40-80	40-80	40-80	40-80	40-80
Salt retention rate [%]	95-98	95-98	95-98	95-98	95-98	95-98
Conductivity [µS/cm]	<20	<20	<20	<20	<20	<20
Storage reservoir [l]	100	100	100	100	100	100
Electrical connection [V]	230 (115)*	230 (115)*	230 (115)*	230 (115)*	230 (115)*	230 (115)*
Power consumption [kW/h]	1.3	1.3	1.3	1.3	1.3	1.3
Power frequency [Hz]	50 (60)*	50 (60)*	50 (60)*	50 (60)*	50 (60)*	50 (60)*
Water intake [diameter]	1/2"	1/2"	1/2"	1/2"	1/2"	1/2"
Concentrate outlet [diameter]	10 mm	10 mm	10 mm	10 mm	10 mm	10 mm
Permeate outflow [diameter, hose]	3/4"	3/4"	3/4"	3/4"	3/4"	3/4"
Height [cm]	81	81	81	81	81	81
Width [cm]	58	58	58	58	58	58
Depth [cm]	58	58	58	58	58	58
Max. water temperature	25°C	25°C	25°C	25°C	25°C	25°C
Max. water pressure	7 bar	7 bar	7 bar	7 bar	7 bar	7 bar
Min. water pressure	2 bar	2 bar	2 bar	2 bar	2 bar	2 bar
Weight (full)	150 kg	151 kg	152 kg	153 kg	154 kg	155 kg
Number of membranes	1	2	3	4	5	6

* RO system 613L--B

** Drinking water quality at 10°C, 3 bar, max. 500 mg/l total salt content ± 5% flow rate

*** Drinking water quality at 500 µS/cm
40% without pre-treatment softening (standard)
80% with pre-treatment softening (option)

1.7 Technical Specifications

Symb.	Description	Type Data
P1	High pressure pump	Procon: 2539/Grundfos MG80B, 12-16.5 bar, 1 × 230 V, 0.55 kW (1 × 115 V, 0.75 kW)*
P2	Transport pump	CHI2-50, 1 × 230 V, 0.80 kW (1 × 115 V, 0.55 kW)*
FI 1	Flow meter	Ø 20 PVC
PI 1	Pressure gauge	0-40 bar, 1/4"
V1	Needle valve (permeate discharge valve)	1/4", Brass
V2	Needle valve (concentrate discharge valve)	1/4", Brass
Y1	Solenoid valve NC	1/2", Brass
PS 1	Pressure controller NO	1/4", 0.5 bar
PS 2	Pressure controller NC	1/2", 2-6 bar
QIS 1 (option)	Conductivity meter	1/2" Connection for sensor

* RO system 613L--B

2. Initial System Operation

2.1 System Placement

System dimensions are 58+58+81 cm (W+D+H). When placing the system, you must ensure that there is enough room to remove the cover. An additional height of 37 cm is required over the system to ensure that the cover can be removed and the system can be rolled out for servicing (e.g. when the system is placed under a table or similar).

The surface upon which it is placed must be able to withstand a load of approx. 160 kg – the possible weight of a full system in operation. The surface must be even and free of grease to allow the RO system to move on its attached rollers.

The system should be set up on a reasonably level surface that is frost-free so that the permeate does not overflow into the supply tank when the container is full.

There should be sufficient space at the back of the system for the water connection, paying special attention to the drain. The drain hose must not be kinked under any circumstances.

There should be sufficient space at the back of the system for the water connection, paying special attention to the drain. The drain hose must not be kinked under any circumstances.

The system should be set up so that the air intake at the front of the system is not covered under any circumstances.

Make sure that the alarm LEDs on the front of the system are visible.

In case of a system malfunction during operation, the level in the tank may rise above maximum. To prevent water damage, the system should be located near a non-pressurized floor drain.

2.2 Water Quality

The raw water to be treated by RO 613L systems must be drinking water quality and may contain a max. of 1000 mg/l TDS.

The raw water may not exceed the following standard values:

Fe:	0.05 mg/l
Mn:	0.02 mg/l
Free chlorine:	0.1 mg/l
Turbidity max.	1.0 NTU
Silt index:	3.0
KMnO ₄ max.:	10 mg/l
Max. temperature:	25°C

The system factory setting for raw water temperature is 10°C.

If you are unsure of the composition of the raw water, a water analysis may be necessary.

The system should be connected to a raw water supply with a pressure of at least 2 bar and a max. of 7 bar.

The treated water (permeate) should have a conductivity of under 20 µmS/cm at 10°C.

2.3 Electrical Connection

All system water connections must meet local regulations and system specifications.

Raw water connection | Supply to the RO system (see section 4.2)

A ball cock must be attached to the input side of the system so that the raw water supply can be interrupted for servicing.

A 1/2" flexible pressure hose should be connected to the system prefilter. The other end of the pressure hose is connected to the raw water supply.

The installation should be made to a supply pipe that is at least 3/4." Using a smaller pipe dimension might increase the risk of system interruptions due to inadequate water pressure and/or a flow rate that is too low, e.g. when flushing membranes during system start-ups.

Connecting the permeate | product water outlet

The supply line to the user (or a subsequent pipe connection) from the transport pump should be made using a 3/4" flexible high pressure hose.

Water that has been completely desalinated can speed up corrosion. For this reason, you should always use corrosion-proof piping for the treated water, e.g. stainless steel or PVC.

Connecting the concentrate outlet hose (see section 4.2)

Remove the "broken" blue hose piece on the concentrate discharge valve and replace with the 10 mm plastic hose supplied.

Please press the hose completely into the coupling.

The end of the hose should lead to the floor drain. The drain hose should not be under the water level in the floor drain; otherwise, there is a risk that the system may draw in the dirty water.

There should not be any kinks or other hindrances to the hoses because this would damage the system membrane(s).

Connecting the overflow hose (see section 4.2)

A 1" plastic hose should be connected from the overflow connection on the back of the system to the floor drain or to another drain below it. This overflow is a safety feature to prevent the system from overflowing if there is a malfunction.

2.4 Electrical Connection

All system electrical connections must meet local legal regulations and system specifications. If these regulations are not followed, this could lead to serious personal injury (electrical shock) and/or permanent system damage.

Should the power connection need to be replaced, please follow the respective wiring diagram (see section 4.3).

Connecting the RO system to the voltage supply (see section 4.2)

The following values apply to the RO 613L system power connection (the values for the RO 613L--B system are in parentheses):

Voltage:	230 Volt–50 Hz (115 Volt–60 Hz)
Fuse:	10 A (32 A)
Highest consumption:	1.5 KW

The system should be connected to 'Phase,' '0' and 'Ground.'

All internal system connections such as the pump control and level control are factory set. You only have to connect the cable from the control to an On|Off switch.

Color of wiring

Blue 0 wire

Brown Phase wire

Yellow | Green Ground wire

2.5 System Startup (see section 4.2)

Before startup, make sure that all water and electrical connections are correct according to our description as well as any local regulations.

Startup procedure

- ▶ Open raw water supply.
- ▶ Make sure that all water connections are leak proof.
- ▶ Remove the short permeate hose from the input to the storage tank and replace it with the supplied hose. Guide the end of the hose to the drain.
- ▶ Open the concentrate discharge valve completely.
- ▶ Turn on the voltage supply of 230 Volt and 50 Hz (115 Volt and 60 Hz for the 613L--B RO system) and activate the main switch on the control.
- ▷ The system is now operational.

During the initial startup, the system must first be rinsed for 20-30 minutes.

After rinsing, the concentrate discharge valve should be adjusted according to the following section.

All system electrical connections must meet local regulations.

2.6 Adjustment of Concentrate Flow Rate

The steps described in this section and the “Setting the Permeate Flow Rate” section must be followed exactly to ensure a correct adjustment of the system.

The correct concentrate flow rate is different for each system depending on the quality of the raw water. Recovery of raw water that is too high can damage the system membranes. Provided that the quality of the raw water is within the prescribed limits, water recovery of 40% can be achieved (with surface water approx. 50%). Water recovery of 70% to a max. of 80% is possible with softened water, depending on the amount of organic material in the water.

The conductivity of the permeate depends on the water recovery of the system, i.e. for a desired salt retention of more than 98%, the water recovery of the system can be set lower. Please note that a salt retention of 98% corresponds to a water recovery of 75%. You can ask your supplier regarding the optimal output amount related to your water quality.

Calculating the System Output Amount (Concentrate)

The output amount is calculated using the following values and formulas:

RO systems	Permeate capacity [l/h]	Concentrate output amount [l/h] (recovery)		
		Ground water	Surface water	Soft water
		40%	50%	75%
613L050	50	75	50	17
613L100	100	150	100	33
613L150	150	225	150	50
613L200	200	300	200	67
613L250	250	375	250	83
613L300	300	450	300	100

The concentrate amount for the systems can be calculated as follows:

$$\text{Concentrate amount [l/h]} = \frac{100 + \text{permeate capacity [l/h]}}{\text{Water recovery (\%)}} - \text{Permeat Leistung [l/h]}$$

Example: RO 613L150 with 40% water recovery

$$\text{Concentrate amount [l/h]} = \frac{100 + 150}{40} - 150 = 225 \text{ l/h}$$

After setting the desired amount of concentrate, fasten the lock nuts of the concentrate valve to ensure that the valve setting is not changed. The concentrate amount should be re-checked after fastening the lock nuts to ensure that the valve has been set properly. Both lock nuts must be attached and secured.

The setting for the concentrate flow rate may not fall below the prescribed value; otherwise, the membranes will be damaged.

2.7 Setting the Permeate Flow Rate

The permeate flow rate is set on the permeate valve. The maximum values for permeate per liter/hour is 50–300 l/h for RO 613L (depending on the RO system type) for a temperature between 10–25°C.

If the temperature should fall below 10°C, reduce the standard flow by approx. 3% per degree Celsius.

For example, if the raw water temperature is 8°C, this means that the permeate flow for the RO 613L200 system will be 6% below the standard flow of 200 l/h, i.e. 188 l/h.

The pressure of pressure gauge must not fall below 12 bar.

In order to reach a standard permeate flow rate, the pressure on the pressure gauge must show between 12 and 16.5 bar.

When the desired pressure and permeate flow rate has been reached, the concentrate amount should be rechecked.

Sartorius Stedim Biotech recommends that both lock nuts be loosened during fine adjustment of the valves.

After both valves have been set, the valves must be secured using the lock nuts. Make sure that you do not change the setting when securing the lock nuts.

After all valves are secured, the system should be turned off and on 4–5 times. Then recheck the flow rate and reset the valves if necessary.

The quality of the permeate|product water should be checked regarding conductivity. If the water quality is <20 µmS/cm, the hose can be attached to the storage tank.

The conductivity meter is available as an accessory.

The 100 l supply tank is now filled with treated water (>20 µmS/cm).

The following operating data should be entered in the included operating log:

- ▶ System operating pressure – read from the high pressure gauge, Operating pressure 12–16.5 bar.
- ▶ Permeate flow rate – read from the flow rate meter, Maximum 50, 100, 150, 200, 250 or 300 l/h.
- ▶ Raw water pressure and raw water temperature should also be entered into the operating log.

3. System Operation

3.1 Automatic

The system control unit is equipped with the following functions:

- Safety shut-down of the transport pump if the storage tank runs dry
- Level sensor to start|stop the booster pump
- Solenoid valve in supply water input
- Pressure switch to start|stop the transport pump (stop 4.0 bar| start 3.7 bar)
- Pump turns off with a time delay of 20 sec.
- Pressure switch for controlling the supply water pressure (> 0.5 bar)
- If the raw water pressure is too low (< 0.5 bar), the system is turned off by the pressure controller after a delay of 15 seconds. The red LED control unit display then indicates a system error.
- Additional level sensor (optional) as overflow protection

3.2 Maintenance (see section 4.2)

The RO system was manufactured and designed with special emphasis on the lowest possible maintenance effort. However, several functions still need to be checked on a regular basis. The maintenance intervals are described in section 3.2.

The following should be checked on a regular basis

The flow rate of treated water should be checked at the flow meter of the RO system.

If the flow rate has fallen more than 10% in relation to the day when the system was set up, see section 3.3 "Troubleshooting."

Check the quality of the treated water using a conductivity meter (option) at the permeate hose during system operation.

If the quality of the treated water exceeds 20 µmS/cm, see section 3.3 "Troubleshooting."

Check the high pressure gauge to determine whether or not the operating pressure is between 12 to 16.5 during system operation. If the pressure is low, see section 3.3 "Troubleshooting."

3.3 Troubleshooting (see section 4.2)

Possible system problems are described in this section.

3.3.1 System Capacity has Dropped

The permeate flow rate also depends on the raw water pressure and the water temperature. If the pressure and temperature are too low, the flow rate is reduced. The flow rate increases when the pressure and temperature increase. Therefore, if the flow rate decrease, it first makes sense to check the pressure and temperature of the raw water before cleaning the membrane.

The flow rate reduction should be read from the flow meter while the system is in operation.

- Check the operating pressure (between 12–16.5 bar) at the pressure gauge during system operation.
- ▷ If operating pressure is low, i.e. < 12 bar, the prefilter should be replaced with a new 5", 5 m prefilter cartridge. The filter should also be replaced if it is very dirty. Also check whether or not the raw water pressure is the same as when the system was first set up. If the pressure has sunk under 2 bar, the error should be searched for in the water supply.

► **Check**

The raw water temperature.

- ▷ If the raw water temperature has decreased compared to the day when the system was set up (winter|summer), the flow rate will also have decreased. At a higher temperature, the flow rate will again increase.
 - ▷ For every °C plus or minus, the flow rate of the system will rise or fall by approx. 3%.
 - ▷ This means that a temperature decrease of 4°C compared to the day the system was set up would mean the flow rate has decreased by around 12%. This is entirely normal and does not require service.
- Rinse the membranes for 20 minutes with the permeate valve completely closed and the concentrate valve completely open.
- ▷ After rinsing is complete, reset the valves as described in section 2.

► **Check**

the drain hose for kinks or any other obstructions.

- ▷ If the flow rate of the system does not increase after correcting this, the membranes are clogged and must be cleaned, see section 3.4 "Cleaning, Replacing Membranes."

► **Check**

the concentrate outlet hose for a kink.

- ▷ Correct the error by replacing the hose.

► **Check**

whether or not the system has been idle for long period of time (1 week or more).

- ▷ Correct the error by running the system for one to two hours and then operate the system a minimum of every three days. Don't forget to reset the concentrate|permeate discharge valve after rinsing is complete, see section 2.

► **Check**

whether or not the seal is intact on the raw water side to make sure that the raw water is not coming into contact with the treated water in the storage tank.

- ▷ Correct any loose seals, empty the storage tank and refill the tank with treated water, <20 µmS/cm.

If you cannot find any of the problems listed here,

- ▷ the problem is with the system membranes, which have to be cleaned and/or replaced, see section 3.4 "Cleaning, Replacing Membranes."

Alarm indicators are displayed via the control box, see also section 3.3.9.

3.3.2 The Quality of the Treated Water is Higher than 20 µmS/cm

3.3.3 Alarm: Lack of Input Pressure

The pressure switch on the supply water input switches the booster pump off when the pressure is too low (< 0.5 bar).

The system then starts 3x automatically before a permanent alarm is generated.

► **Check**

whether or not the prefilter is clogged.

- ▷ Correct the error by turning off the raw water supply and reducing the pressure on the prefilter. Then replace the prefilter with a new 5", 5 µm prefilter cartridge. Turn off the power to the system for 20 sec. and then turn it on again.

- ▶ **Check**
the raw water for the required pressure.
- ▶ **Look for the error in the raw water supply.**
- ▷ When the raw water pressure has been restored, the system power supply should be turned off and then back on again after 20 sec.

If this does not correct the problem,

- ▷ it is possible that the pressure controller at the RO system prefilter is defective or the PC board in the control unit is defective.

3.3.4 Alarm: Transport Pump

- ▶ **Check**
that water is being fed by the transport pump.
- ▶ **Connect a water user to the outlet of the transport pump.**
- ▷ If the pump operates, then the cause of the problem is somewhere else.
The system and the pump control are OK.

3.3.5 Alarm: Tank Fill Level Too Low

- ▶ **Check**
whether or not the tank has run dry.
- ▶ **Let the tank fill up completely, then the transport pump will turn on automatically.**

When the tank empties completely, the transport pump will only turn on when the container is completely full. This is controlled by the system level regulator.

- ▶ **Check**
whether or not the transport pump pressure switch is damaged.
- ▶ **Bypass the pressure switch by bridging both plugs.**
- ▷ If the pump runs after this and water is supplied, the cause is the pressure switch, which must be replaced.
- ▶ **Check**
whether or not the function switch of the transport pump is set to Off.
- ▷ If none of these problems are the cause, then the control board or the transport pump is damaged and replacement is necessary.

3.3.6 Alarm: Tank Fill Level Too High

- ▶ **Check**
whether or not the solenoid valve is closed in the input.
- ▷ If it continues to drip from the permeate | product water hose when the system is inactive, then the solenoid valve is defective and must be replaced.
- ▶ **Check**
whether or not the level sensor is blocked.
- ▷ If the level sensor is not blocked or the blockage cannot be removed, please contact your Sartorius Stedim Biotech service technician.

3.3.7 System Will Not Start

- ▶ **Check**
whether or not the main power supply is connected.
- ▶ **Check**
whether or not the controller main switch is turned on.
- ▶ **Check**
whether or not the system should operate at all.
→ Full storage tank or no "requirement" for water.
- ▷ If none of the above problems are the cause, then there is a problem with the high pressure pump or the control unit. These should be checked.


3.3.8 Transport Pump Turns On | Off

- The system transport pump turns on and off every 10 to 15 seconds without treated water being used.
- **Check**
whether or not there is a leak in the connection between the system output and the point of consumption of the treated water or that there is damage at the point of consumption.
 - ▷ A damaged or leaky valve could trigger water consumption when none is required, which would cause the reservoir pump to constantly turn on and off.
 - **Check**
whether or not the return valve on the intake rod of the transport pump is leaking and | or damaged.
 - ▷ The valve should be replaced if leaking and | or damaged.
 - **Check**
whether or not the pressure tank for the transport pump contains enough air.
 - ▷ The pressure tank must have a preliminary pressure of 2.9 bar. If not, please contact your Sartorius Stedim Biotech service technician.

3.3.9 Alarm Display

Alarm functions:

Booster pump stops: (high pressure pump operating normally)	(1) lights + (Power)	(3) flashes slowly (High Pressure Pump)
Transport pump stops: (booster pump not operating)	(1) lights + (Power)	(4) flashes slowly (Transport Pump)
Transport pump ON alarm: (system not operating at all - continual status)	(1) lights + (Power)	(6) lights + (4) flashes slowly (Alarm) (Transport Pump)
Transport pump restart: (system not operating at all - continual status)	(1) lights + (Power)	(6) lights + (4) flashes quickly (Alarm) (Transport Pump)
Input pressure too low (restart): (system not operating at all - continual status)	(1) lights + (Power)	(2) flashes slowly (Inlet)
Input pressure too low (alarm): (system not operating at all - continual status)	(1) lights + (Power)	(6) lights + (2) flashes quickly (Alarm) (Inlet)
Tank fill volume too low (alarm): (transport pump stops - high pressure pump operating normally)	(1) lights + (Power)	(5) flashes slowly (Level)
Tank fill volume too high (alarm): (system not operating at all - continual status)	(1) lights + (Power)	(6) lights + (5) flashes quickly (Alarm) (Level)


3.3.10 Time Setting

It is possible to manually adjust the various pump time settings such as start, stop, alarm and start delay.

Bridge|Jumper function:

1. Booster pump
2. Transport pump
3. RO module selection
- 4.|5.|6. Transport pump

Factory settings:

Delay Off 5 seconds		Delay On 15 seconds
ALARM Off No alarm		ALARM On Stops after 20 minutes
MODEL Off RO		MODEL On RO
		2 On 30 seconds starts max. 60/h
		4 On 30 seconds starts max. 90/h
		6 On 30 seconds starts max. 109/h
		16
		32

3.4 Cleaning|Replacing Membranes

3.4.1 Cleaning Procedure

Do the following before cleaning:

- ▶ Remove approx. 50 liters of treated water from the storage tank.
- ▶ Disconnect the system power supply.
- ▷ Disconnect the raw water supply to the system.

Cleaning should be carried out as follows:

- ▶ Remove the water inlet hose from the raw water intake.
- ▶ Open the prefilter housing and remove the prefilter cartridge.
Fill the filter housing with raw water.
- ▶ Open the concentrate discharge valve completely and close the permeate discharge valve.
- ▶ Take a plastic bucket (min. 25 l) and fill it with 25 l of warm raw water (approx. 40°C).
- ▶ The container must be higher than the prefilter housing.
- ▶ Remove the permeate hose from the storage tank and put it in the container with the cleaning liquid (see below).
- ▶ Put the raw water supply hose into the container all the way to the bottom.
- ▶ Bypass the wires of the attached pressure switch (behind the prefilter).
Do this by crossing the 2 wires or bridging the 1 and 4 terminals.
- ▶ Put the concentrate hose in the container with the warm water.
- ▶ Start the system by connecting the power. The system will start operation at a somewhat higher noise level until the warm water has been added.

Operate the system in this way approx. 20–30 minutes until the system components (membranes | pumps) have reached a temperature of approx. 40°C. The water in the container should be replaced on a continual basis to keep the temperature at approx. 40°C during the warming period.

- ▶ When the system components have reached a temperature of approx. 40°C, turn off the system power supply.
- ▶ Re-fill the container with warm water (40°C). We recommend using treated water for this purpose.
- ▶ Add 2% by volume of citric acid to the 25 liters of water, i.e. 0.5 kg citric acid to 25 liters of water.
- ▶ Restart the system by connecting the power.
- ▶ Operate the system for 5 minutes with this citric acid solution.
- ▶ Stop the system for 5 minutes.
- ▶ Then restart the system for 5 minutes. Repeat this process 3 times, i.e. 3 × 5 minutes in operation, after each time in operation resting for 5 minutes.
- ▶ If the membranes are really clogged, an additional cleaning cycle may be necessary with a new citric acid solution. This additional cleaning should not be required if the capacity has only been reduced by 10%.
- ▶ Once cleaning is complete, the cleaning liquid should be disposed off in a proper manner. Then reconnect the raw water hose to the raw water supply.
- ▶ Reattach the concentrate hose to the outlet.
- ▶ Reattach the permeate hose to the outlet.
- ▶ Insert the 5", 5 µm prefilter cartridge (possibly a new one) into the prefilter housing and reconnect the wires to the pressure switch.
- ▶ Flush the system approx. 20–30 minutes.
- ▶ Once the flushing is complete, reset the discharge valve | recirculation valve, see section 6 System Startup.

- ▶ When the system pressure gauge has reached 12.0–16.5 bar, check at the permeate hose whether or not the water quality is $<20\text{ }\mu\text{mS/cm}$ (a conductivity meter is available as an option). If the water quality is $>20\text{ }\mu\text{mS/cm}$, operate the system for 10–20 minutes before rechecking.
- ▶ Check whether or not the flow rate at the flow meter is satisfactory.
- ▶ When the flow rate and water quality of the system are satisfactory, reattach the permeate hose to the tank.
- ▷ The system is now ready for normal operation.

If the system does not achieve full capacity after cleaning, you should run through the cleaning procedure once more. If after this additional cleaning, the performance of the membrane(s) is still unsatisfactory, the membrane(s) should be replaced.

3.4.2 Tank Cleaning Procedure

Do the following before cleaning:

- ▶ Remove the permeate/product water until the removal activates the RO high pressure pump.
- ▶ Stop removing the water and let the system run until the maximum volume is reached. The high pressure pump turns off automatically.
- ▷ Then turn off the system.

Cleaning should be carried out as follows:

- ▶ Add H_2O_2 to the product water until you have created a 1% H_2O_2 cleaning solution.
 - ▶ Let the cleaning solution have a chance to work for 1 hour.
 - ▶ After the application time remove the permeate hose from the tank. Put this and the product water discharge hose into the non-pressurized drain.
 - ▶ Turn the system on and empty the tank completely (step "A").
 - ▶ Then turn the system off.
 - ▶ Reattach the permeate hose to the tank (step "C").
 - ▶ Turn the system on and let the tank fill completely (step "D").
 - ▶ Turn the system off.
- Remove the permeate hose after the tank application time.
Put this and the product water discharge hose into the non-pressurized drain (step "E").
- ▷ Repeat this process (steps "A" through "E" three times and complete the third time with step "D").

3.4.3 Replacing Membranes

The steps should be carried out before replacing the membranes:

- ▶ Remove approx. 50 liters of treated water from the storage tank.
- ▷ Disconnect the system from the power supply.

Membrane replacement should be carried out as follows:

Note the respective positions of the hoses before removal to ensure proper reassembly.

- ▶ Remove the plastic hoses at the end of the membrane tubes.
- ▶ The hoses can be removed by pressing down on the orange-colored ring.
- ▶ Remove the safety pin and then remove the U fastener at the end of the membrane tube (the U fastener secures the position of the end pieces).

- ▶ Carefully remove the end pieces up out of the membrane tube by using a sideways back and forth motion while pulling at the same time.

Note the end of the membrane(s) with the large black lip seal. When installing the new membrane(s), this seal must be in the same position (top or bottom).

- ▶ Remove the membrane(s) from the membrane tube.
- ▶ After membrane replacement, reinsert the end pieces and secure with the U fastener and pin. Then reattach all hoses.

The “ring” should be secured tightly when sliding in the hose connection and the hose should be inserted completely.

- ▶ After all connections are complete and the end pieces secured with the U fastener, the system should be reconnected to the power supply.
- ▶ The prefilter may need to be replaced with a new 5", 5 µm prefilter cartridge.
- ▶ Reconnect the raw water.
- ▶ Open the discharge valve completely.
- ▶ Close the recirculation valve completely.
- ▶ Remove the permeate hose and run it to the drain.
- ▶ Reconnect the system power supply.
- ▶ Flush the system in this way for 20 to 30 minutes.
- ▶ Set the concentrate discharge valve and permeate discharge valve.
- ▶ Check the operating pressure on the pressure gauge to determine whether or not normal operating pressure has been achieved (12–16.5 bar).
- ▶ Check that the water quality is within specifications (standard 20 µmS/cm). This can be read from the conductivity meter (optional). If the water quality is satisfactory, connect the hose to the tank.
- ▶ Check the flow meter for the desired system flow rate.
- ▷ The system is now operational and can be operated normally.


The following information should be entered in the operating log:

- Date of membrane replacement
- New system flow rate (flow meter)
- Water quality (µmS/cm)
- System operating pressure (pressure gauge)
- Raw water temperature
- Raw water pressure
- Tank cleaning


For your notes:

4. Appendix


4.1 Flow Diagram

4.2 Detailed System Drawings:
4.2.1 Layout Overview


4.2.2 Front Side View

Key


- A: Air intake for electrical components
- B: Handle for system transport
- C: Wheels
- D: Flow meter for reading permeate flow rate
- E: LED alarm

4.2.3 Back Side View

- Key**
- A: Supply water input (1/2" pressure hose)
 - B: Permeate|Product water outlet (3/4" pressure hose)
 - C: Overflow connection (must be fed to the floor drain)


4.2.4 Top View


Key

- | | |
|--------------------------------|---------------------|
| A: Permeate hose | H: Level sensor rod |
| B: Concentrate discharge valve | I: Transport pump |
| C: Control box | J: RO membrane(s) |
| D: Permeate discharge valve | K: Pressure switch |
| E: Pressure gauge | L: Pressure tank |
| F: Flow meter | M: Prefilter |
| G: High pressure pump | |

4.3 Wiring Diagram


Power

¹ Capacitor: 20 µf / 450 VDB

4.4 Main Board

- POWER
green LED
- Input valve Y1
green LED
- RO pump P1
green LED
- Transport pump P2
green LED
- Product water level
yellow LED
- Error
red LED

For your notes:

4.5 Service and Maintenance of the RO System

Service and Maintenance of the RO System	Every day	Every week	Every 5-8 weeks	Semi-annually (or as needed)	Every 52nd week
(not standard) Check soft water supply (<0.5°dH).	×				
Check permeate flow rate at flow meter (see page 18, F) while system is in operation.		×			
Check UV lamp*.					
Check system operating pressure. On pressure gauge (see page 18, E) while system is in operation.			×		
Check and clean the system and the tank.				×	
Clean membrane(s). When flow rate is reduced by 10%.					×

For your notes:

4.6 Operating Log

[illegible]

For your notes:

4.7 Disposal

The packaging is to be taken to a local waste disposal site if no longer required. The packaging comprises of environmentally-friendly materials that can be used as secondary raw materials.

The device, including accessories and batteries, is not to be thrown into the household waste. EU legislation in Member States requires electrical and electronic equipment to be collected separately from unsorted municipal waste so that it may be recycled.

In Germany and several other countries, Sartorius Stedim Biotech GmbH itself assumes responsibility for the return and conformant disposal of its electronic and electrical products. These products may not be placed with household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other member nations of the European Economic Area (EEA), please contact our local service technicians or our Service Center in Goettingen, Germany:

Sartorius Stedim Biotech GmbH
Biotechnology Service Center
August-Spindler-Strasse 11
37079 Goettingen, Germany

In countries that are not members of the European Economic Area (EEA) or where no Sartorius Stedim Biotech subsidiaries or dealerships are located, please contact your local authorities or a commercial disposal operator.

Remove the batteries and hand them in to a collection point prior to disposal | scrapping of the device.

Sartorius AG, its affiliates, subsidiaries, dealers and distributors will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the package leaflet/the included CD ROM or to our web site (www.sartorius-stedim.com) for more detailed information regarding addresses for repair service or disposal of your device.

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U.K.

Sartorius Stedim UK Limited
Longmead Business Park
Blenheim Road, Epsom
Surrey KT19 9 QQ
Phone +44.1372.737159
Fax +44.1372.726171

America

USA

Sartorius Stedim North America Inc.
5 Orville Drive
Bohemia, NY 11716
Toll-Free +1.800.368.7178
Fax +1.631.254.425

Sartorius Stedim SUS Inc.
1910 Mark Court
Concord, CA 94520
Phone +1.925.689.6650
Toll Free +1.800.914.6644
Fax +1.925.689.6988

Sartorius Stedim Systems Inc.
201 South Ingram Mill Road
Springfield, MO 65802
Phone +1.417.873.9636
Fax +1.417.873.9275

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Asia|Pacific

China

Sartorius Stedim Beijing
Representative Office
No. 33, Yu'an Road,
Airport Industrial Zone B,
Shunyi District
Beijing 101300
Phone +86.10.80426516
Fax +86.10.80426580

Sartorius Stedim Shanghai
Representative Office
Room 618, Tower 1, German Centre,
Shanghai, PRC., 201203
Phone +86.21.28986393
Fax +86.21.28986392.11

Sartorius Stedim Guangzhou Office
Room 704, Broadway Plaza,
No. 233-234 Dong Feng West Road
Guangzhou 510180
Phone +86.20.8351.7921
Fax +86.20.8351.7931

India

Sartorius Stedim India Pvt. Ltd.
10, 6th Main, 3rd Phase Peenya
KIADB Industrial Area
Bangalore – 560 058
Phone +91.80.2839.1963|0461
Fax +91.80.2839.8262

Japan

Sartorius Stedim Japan K.K.
KY Building, 8-11
Kita Shinagawa 1-chome
Shinagawa-ku
Tokyo 140-0001
Phone +81.3.3740.5407
Fax +81.3.3740.5406

Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur
Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore

Sartorius Stedim Singapore Pte. Ltd.
10, Science Park Road, The Alpha
#02-25, Singapore Science Park 2
Singapore 117684
Phone +65.6872.3966
Fax +65.6778.2494

Australia

Sartorius Stedim Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175
Phone +61.3.8762.1800
Fax +61.3.8762.1828

Sartorius Stedim Biotech GmbH
August-Spindler-Str. 11
37079 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part of
this publication may be reprinted
or translated in any form or by any
means without the prior
written permission of Sartorius
Stedim Biotech GmbH.
The status of the information,
specifications and illustrations
in this manual is indicated by
the date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make changes
to the technology, features,
specifications and design of the
equipment without notice.

Status:
July 2009,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany